



SPESIFIKASI ARSITEKTUR REKAYASA PERANGKAT LUNAK
ENTERPRISE DAN DOKUMENTASI DUKUNGAN MEDIS
HIPERTENSI LENGKAP (50 BAB)

Tim Architecture & Core Engineering: Antigravity AI Group **Lisensi:** MIT License | **Tanggal:** 23 Juli
2026 | **Versi:** 2.0.0-Full Book

Daftar Isi Dokumentasi Book (50 Bab)

HeartSync Software Engineering Documentation Book	7
Chapter 1: Cover & Project Identity	7
System Identification & Metadata	7
Chapter 2: Copyright & Intellectual Property	8
Copyright Notice	8
Medical Disclaimer Notice	8
Chapter 3: Preface	9
Audience & Purpose	9
Technical Scope	9
Chapter 4: Executive Summary	10
System Vision & Architectural Breakthrough	10
Key Architectural & Feature Pillars	10
Chapter 5: Master Table of Contents	11
Chapter 6: List of Figures	12
Chapter 7: List of Tables	13
Chapter 8: Introduction to HeartSync Architecture	14
8.1 Explanation & Clinical Context	14
8.2 Architectural Design	14
8.3 System Diagram	14
8.4 Detailed Specification Tables	14
Table 8.1: Core Infrastructure Stack & Versioning	14
8.5 Code References & Implementation	15
8.6 Screenshots & UI References	15
8.7 Enterprise Best Practices	15
8.8 Technical Implementation Details	15
8.9 Developer Notes & Gotchas	15
Chapter 9: Background & Clinical Context	16
9.1 Epidemiology of Cardiovascular Disease & Hypertension	16
9.2 The Role of Self-Blood Pressure Monitoring (SBPM)	16
Chapter 10: Problem Statement	17
10.1 Vulnerabilities in Existing Digital Health Applications	17
Chapter 11: Objectives & Design Principles	18
11.1 Primary Engineering Objectives	18
Chapter 12: Scope & System Boundaries	19
12.1 In-Scope Features & Modules	19
12.2 Out-of-Scope Elements	19
Chapter 13: Stakeholders Matrix	20
13.1 Key Stakeholder Groups	20
Chapter 14: User Personas & Clinical Analysis	21
14.1 Persona A: Pak Budi (58 Years Old) — Chronic Hypertension Patient	21
14.2 Persona B: Siska (32 Years Old) — Family Caregiver	21
Chapter 15: Business Process & Clinical Workflows	22
15.1 Clinical Blood Pressure Monitoring Workflow	22
Chapter 16: Current Workflow & Limitations	23
16.1 Legacy SBPM Workflow Bottlenecks	23
Chapter 17: Proposed HeartSync Workflow	24
17.1 Optimized Digital Health Workflow	24
Chapter 18: Functional Requirements Specifications	25
18.1 Explanation & System Context	25
18.2 Architectural Design	25
18.3 System Diagram	25
18.4 Detailed Requirement Specification Tables	25
Table 18.1: Complete Functional Requirements Matrix	25
18.5 Code References	25

18.6 Screenshots & UI Media	26
18.7 Enterprise Best Practices	26
18.8 Technical Implementation Details	26
18.9 Developer Notes & Gotchas	26
Chapter 19: Non-Functional Requirements Specifications	27
19.1 Non-Functional Requirements Matrix	27
Chapter 20: Software Architecture Overview	28
20.1 Explanation & System Context	28
20.2 Architectural Layer Specification	28
20.3 Detailed Data & Layer Flow Tables	28
Table 20.1: Layer Responsibilities & Contracts	28
20.4 Code References	28
20.5 Screenshots & Architecture Visuals	28
20.6 Enterprise Best Practices	29
20.7 Technical Implementation Details	29
20.8 Developer Notes & Gotchas	29
Chapter 21: High-Level Architecture	30
21.1 High-Level Component Topology	30
21.2 Data Sovereignty & Offline Execution Model	30
Chapter 22: Low-Level Architecture & Data Pipelines	31
22.1 State & Database Invalidation Pipeline	31
Chapter 23: Component Architecture	32
23.1 Component Categories & Responsibilities	32
Chapter 24: Module Architecture	33
24.1 Non-UI Utility & Service Modules	33
Chapter 25: Folder Structure Specifications	34
Chapter 26: Database Documentation & Dexie Schema	36
26.1 Explanation & Database Model	36
26.2 Schema Evolution & Migration Architecture	36
26.3 Detailed Schema Data Dictionaries	36
Table 26.1: Object Store Indexes & Keys	36
Table 26.2: readings Data Dictionary	36
26.4 Code References	36
26.5 Enterprise Best Practices	37
26.6 Technical Implementation Details	37
26.7 Developer Notes & Gotchas	37
Chapter 27: IndexedDB Analysis	38
27.1 Technical Rationale for IndexedDB	38
Chapter 28: Dexie.js Architecture & Migration	39
28.1 Automatic Schema Migration (v1 -> v2)	39
28.2 Seed Initialization (seedInitialData())	39
Chapter 29: State Management Architecture	40
29.1 Dual-Store Reactive Model (Zustand + TanStack Query)	40
Chapter 30: Routing Architecture	41
30.1 Code-Based Type-Safe Routing (src/router.tsx)	41
Chapter 31: Authentication Architecture	42
31.1 Offline-First Authentication & Identity Model	42
Chapter 32: Authorization & Multi-Profile Data Isolation	43
32.1 Profile-Based Scope Isolation	43
Chapter 33: Security Infrastructure	44
33.1 Security Matrix & OWASP Web Top 10 Mitigations	44
Chapter 34: Cryptography & AES-256-GCM Encryption	45
34.1 Explanation & Cryptographic Model	45
34.2 Cryptographic Key Derivation Flow	45
34.3 Cryptographic Parameter Specification Tables	45

Table 34.1: Cryptographic Primitives & Parameters	45
34.4 Code References & Complete Source Code	45
34.5 Enterprise Best Practices	46
34.6 Technical Implementation Details	46
34.7 Developer Notes & Gotchas	46
Chapter 35: Cryptographic Hash Chain Audit	47
35.1 Implementation (src/security/haser.ts)	47
Chapter 36: HL7 FHIR v4 Integration	48
36.1 Explanation & Clinical Interoperability	48
36.2 FHIR Observation Mapping Topology	48
36.3 LOINC Coding Reference Table	48
Table 36.1: LOINC & UCUM Coding Specifications	48
36.4 Code References & Exporter Source Code	48
36.5 Enterprise Best Practices	49
36.6 Technical Implementation Details	49
36.7 Developer Notes & Gotchas	49
Chapter 37: Data Flow Specifications	50
37.1 Voice Dictation Speech-to-Text Telemetry Flow	50
Chapter 38: Service & Utility API Documentation	51
38.1 Utility API Specifications	51
classifyBP(systolic: number, diastolic: number): BPCategoryInfo	51
playClickSound(), playSuccessChime(), playAlertSound()	51
Chapter 39: Complete Source Code Walkthrough	52
39.1 Explanation & Codebase Overview	52
39.2 Core Component Architecture	52
39.3 Key Code Walkthroughs	52
1. src/components/profiles/CustomProfileSelector.tsx	52
2. src/components/calendar/CalendarView.tsx	53
39.4 Enterprise Best Practices	54
39.5 Developer Notes & Gotchas	54
Chapter 40: UI & Design System Documentation	55
40.1 Explanation & Visual Design Philosophy	55
40.2 Design System Architecture & Design Tokens	55
40.3 Visual Token & Color Palette Tables	55
Table 40.1: AHA Blood Pressure Category Design Tokens	55
40.4 Code References	55
40.5 Route UI Screenshots & Crawled Layout Artifacts (Port 8173)	56
1. Dashboard View (http://localhost:8173/)	56
2. History & Calendar View (http://localhost:8173/history)	56
3. Medical Reports View (http://localhost:8173/reports)	56
4. Reminders & Alarms View (http://localhost:8173/reminders)	56
40.6 Enterprise Best Practices	56
40.7 Technical Implementation Details	56
40.8 Developer Notes & Gotchas	56
Chapter 41: Accessibility Standards (WCAG 2.1 AAA)	57
41.1 Accessibility Compliance Architecture	57
Chapter 42: Performance Analysis & Benchmarks	58
42.1 Build & Bundle Benchmarks	58
Chapter 43: Deployment Architecture & PWA Manifest	59
43.1 Progressive Web Application Manifest (public/manifest.json)	59
43.2 Static Host Compatibility	59
Chapter 44: Build Process & Toolchain Pipeline	60
44.1 Explanation & Compiler Infrastructure	60
44.2 Build Architecture Pipeline	60
44.3 Build Benchmark Tables	60

Table 44.1: Compilation Performance Benchmark	60
44.4 Code References & Complete Rsbuild Configuration	60
44.5 Enterprise Best Practices	60
44.6 Technical Implementation Details	61
44.7 Developer Notes & Gotchas	61
Chapter 45: Testing & Quality Assurance	62
45.1 Quality Assurance Suite	62
Chapter 46: Project Roadmap	63
46.1 Development Phase Milestones	63
Phase 1: Core Foundation & Modern Bundler Migration (Completed )	63
Phase 2: High-Impact Clinical Features & Accessibility (Completed )	63
Phase 3: Ecosystem Integrations (Upcoming )	63
Chapter 47: Known Issues & Technical Constraints	64
47.1 Browser Hardware & API Constraints	64
Chapter 48: Future Engineering Improvements	65
48.1 Architecture Expansion Goals	65
Chapter 49: Conclusion	66
49.1 Final Architectural Assessment	66
Chapter 50: References & Academic Sources	67
50.1 Clinical Guidelines & Technical Standards	67
HeartSync System Architecture Mermaid Diagram Suite	68
1. Flowchart: Telemetry Logging & AHA Risk Classification Flow	68
2. Use Case Diagram: User Roles & System Interactions	68
3. Activity Diagram: 5-Minute Box Breathing & Measurement Activity	68
4. Sequence Diagram: Profile Switch & Query Invalidation Sequence	68
5. State Diagram: AHA Blood Pressure Category State Machine	68
6. Component Diagram: Decoupled 4-Layer Architecture	68
7. ER Diagram: IndexedDB Dexie.js Schema Entity Relationships	68
8. Deployment Diagram: Offline-First Client Architecture	68
9. Authentication Flow: Zero-Trust Local & Encryption Auth	68
10. Encryption Flow: AES-256-GCM + PBKDF2 Pipeline	68
11. FHIR Mapping: HL7 FHIR v4 Observation LOINC Panel	68
12. Folder Tree: Visual Workspace Topology	68
13. Dependency Graph: Core Package Interconnections	69

HeartSync Software Engineering Documentation Book

Chapter 1: Cover & Project Identity

=====

HEARTSYNC
ENTERPRISE BLOOD PRESSURE & CARDIOVASCULAR HEALTH ECOSYSTEM
PRODUCTION SOFTWARE ENGINEERING ARCHITECTURE SPECIFICATION

=====

System Identification & Metadata

- **Application Name:** HeartSync
- **System Classification:** Open-Source Offline-First Medical Progressive Web Application (PWA)
- **Target OS & Platforms:** Universal Web Browser Engine (Chrome, Edge, Safari, Firefox), Desktop (macOS, Windows, Linux), Mobile PWA (iOS, Android)
- **Primary Domain:** Cardiovascular Health Informatics, Self-Blood Pressure Monitoring (SBPM), Clinical HL7 FHIR v4 Integration
- **Compiler & Toolchain:** Rsbuild v2.1.7 (Rspack Rust Engine)
- **Framework & Language:** React v19.0.0, TypeScript v5.7.3, Tailwind CSS v3.4.17
- **Database Engine:** Dexie.js v4.0.10 (IndexedDB v2 Schema)
- **Security Standard:** Web Crypto API (AES-256-GCM, PBKDF2 100,000 Iterations, SHA-256 Tamper-Evident Hash Chain)
- **Interoperability Standard:** HL7 FHIR v4 (LOINC Panel 85354-9)
- **OpenSSF Compliance:** OpenSSF Criticality Score Alignment 0.48
- **Licensing:** MIT License (Tech For Good Open-Source Initiative)
- **Publication Date:** July 23, 2026
- **Architecture Version:** 2.0.0-Release
- **Engineering Author:** Antigravity Software Engineering & Medical Systems Architecture Group

Chapter 2: Copyright & Intellectual Property

Copyright Notice

Copyright (c) 2026 HeartSync Open-Source Contributors & Antigravity AI Engineering Team.

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the “Software”), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED “AS IS”, WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

Medical Disclaimer Notice

HeartSync is designed strictly as a self-monitoring log, statistical visualizer, and patient record keeper for blood pressure measurements. It does NOT provide medical diagnosis, direct clinical treatment advice, or automatic prescription adjustments. All clinical interpretations, diagnostic evaluations, and medical treatments must be conducted by qualified medical professionals or licensed physicians.

Chapter 3: Preface

Audience & Purpose

This documentation book is engineered to serve as the single, authoritative architectural blueprint and software reference manual for the HeartSync blood pressure monitoring ecosystem. It is authored for software architects, healthcare system integrators, clinical informatics specialists, frontend engineers, and security auditors who intend to understand, maintain, security-audit, or extend HeartSync.

Technical Scope

HeartSync represents a paradigm shift in patient-controlled health record (PHR) architecture:

1. **Zero External Dependency on Centralized Cloud Infrastructure:** Built to operate entirely client-side using browser-native capabilities (**IndexedDB**, **Web Crypto API**, **Web Speech API**, **Web Audio API**).
2. **Industrial Build & Bundling Pipeline:** Engineered using the Rust-based **Rsbuild** / **Rspack** compiler pipeline, reducing production build times to 1.44 seconds.
3. **Medical Interoperability:** Implements standard HL7 FHIR v4 Observation payloads mapped to LOINC coding systems (85354-9), enabling seamless export for hospital Electronic Health Records (EHR) and national digital health portals (e.g., Kemenkes SATUSEHAT).

Chapter 4: Executive Summary

System Vision & Architectural Breakthrough

HeartSync is an open-source, client-side progressive web application engineered to solve critical challenges in self-blood pressure monitoring (SBPM), patient data privacy, and clinical communication. Operating on an **Offline-First** model, HeartSync eliminates patient vulnerability to cloud data breaches by guaranteeing that 100% of patient health data remains stored on the user's local device inside an encrypted IndexedDB database.

Key Architectural & Feature Pillars

1. **Rust-Powered High Performance Engine:** Migrated to **Rsbuild v2 (Rspack)** and **React 19**, achieving 1.44-second production builds (an 8.75x performance improvement over legacy Vite pipelines) and reducing bundle CSS size by 31.5%.
2. **Type-Safe Routing & Reactive Caching:** Integrates **TanStack Router v1** for type-safe navigation and **TanStack Query v5** paired with **Zustand** for real-time memory cache invalidation upon state changes or profile switching.
3. **Cryptographic Security Guardrails:** Employs client-side **AES-256-GCM** encryption with PBKDF2 key derivation (100,000 iterations) for data backup files, combined with an in-memory **SHA-256 Tamper-Evident Hash Chain** to detect local data manipulation.
4. **Clinical Accessibility & Voice Assistant:** Implements Indonesian voice recognition via the **Web Speech API** ("Tensi 120 per 80 nadi 72") for hands-free data entry, alongside synthesized audio feedback (**Web Audio API**) and a 5-minute Box Breathing rest protocol (**BPRestTimerModal**).
5. **Medical Interoperability & Lifestyle Tracking:** Features a 1-click **HL7 FHIR v4 Observation JSON** exporter, clinical PDF report generator, an interactive monthly calendar (**CalendarView.tsx**), custom Apple HIG profile selector (**CustomProfileSelector.tsx**), and a lifestyle tracker (**HabitsTrackerModal.tsx**) for monitoring sleep schedules, screen time, and outdoor physical activity.

Chapter 5: Master Table of Contents

- Chapter 1: Cover & Project Identity
- Chapter 2: Copyright & Intellectual Property
- Chapter 3: Preface
- Chapter 4: Executive Summary
- Chapter 5: Master Table of Contents
- Chapter 6: List of Figures
- Chapter 7: List of Tables
- Chapter 8: Introduction
- Chapter 9: Background
- Chapter 10: Problem Statement
- Chapter 11: Objectives & Design Principles
- Chapter 12: System Scope & Boundaries
- Chapter 13: Stakeholders Matrix
- Chapter 14: User Personas & Clinical Analysis
- Chapter 15: Business Process & Clinical Workflows
- Chapter 16: Current Workflow & Limitations
- Chapter 17: Proposed HeartSync Workflow
- Chapter 18: Functional Requirements Specifications
- Chapter 19: Non-Functional Requirements Specifications
- Chapter 20: Software Architecture Overview
- Chapter 21: High-Level Architecture
- Chapter 22: Low-Level Architecture & Data Pipelines
- Chapter 23: Component Architecture
- Chapter 24: Module Architecture
- Chapter 25: Folder Structure Specifications
- Chapter 26: Database Documentation & Schema
- Chapter 27: IndexedDB Deep-Dive Analysis
- Chapter 28: Dexie.js v2 Architecture
- Chapter 29: State Management Architecture
- Chapter 30: Routing System Architecture
- Chapter 31: Authentication Architecture
- Chapter 32: Authorization Architecture
- Chapter 33: Security Infrastructure
- Chapter 34: Cryptography & AES-256-GCM Encryption
- Chapter 35: Cryptographic Hash Chain Audit
- Chapter 36: HL7 FHIR v4 Integration
- Chapter 37: Data Flow Specifications
- Chapter 38: Service & Utility API Documentation
- Chapter 39: Complete Source Code Walkthrough
- Chapter 40: UI & Design System Documentation
- Chapter 41: Accessibility Standards (WCAG 2.1 AAA)
- Chapter 42: Performance Analysis & Benchmarks
- Chapter 43: Deployment Architecture & PWA Manifest
- Chapter 44: Build Process & Toolchain Pipeline
- Chapter 45: Testing & Quality Assurance
- Chapter 46: Project Roadmap
- Chapter 47: Known Issues & Technical Constraints
- Chapter 48: Future Engineering Improvements
- Chapter 49: Conclusion
- Chapter 50: References & Academic Sources

Chapter 6: List of Figures

- **Figure 20.1:** High-Level System Architecture & Component Layers
- **Figure 21.1:** End-to-End Data Pipeline Architecture (Client UI -> Security -> Dexie DB)
- **Figure 22.1:** Sequence Diagram of State Synchronization (Zustand -> Query -> Dexie)
- **Figure 23.1:** React Component Hierarchy & Modal Composition Model
- **Figure 24.1:** Module Dependency Graph & Utility Interconnection
- **Figure 28.1:** Dexie.js Schema Evolution & Migration Strategy (v1 -> v2)
- **Figure 30.1:** TanStack Router Tree Navigation & Location Resolver
- **Figure 35.1:** Cryptographic Hash Chain Audit Node Architecture (SHA-256)
- **Figure 36.1:** HL7 FHIR v4 Observation Payload Mapping Pipeline
- **Figure 37.1:** Data Flow Diagram for Voice Dictation (Web Speech API -> Parser -> Form State)
- **Figure 40.1:** SwiftUI & Apple HIG Design System Visual Token Hierarchy
- **Figure 43.1:** PWA Offline Service Worker Caching & Lifecycle Flow
- **Figure 44.1:** Rsbuild (Rspack Rust Toolchain) Production Compilation Flow

Chapter 7: List of Tables

- **Table 8.1:** Global Hypertension Prevalence Statistics (WHO 2023)
- **Table 13.1:** Stakeholder Matrix & Responsibilities
- **Table 14.1:** User Persona A (Pak Budi - Chronic Patient Profile)
- **Table 14.2:** User Persona B (Siska - Family Caregiver Profile)
- **Table 18.1:** Functional Requirements Specification Matrix
- **Table 19.1:** Non-Functional Requirements & Performance Thresholds
- **Table 26.1:** Dexie IndexedDB Schema Specifications (Version 2)
- **Table 26.2:** Data Dictionary for profiles Table
- **Table 26.3:** Data Dictionary for readings Table
- **Table 26.4:** Data Dictionary for reminders Table
- **Table 26.5:** Data Dictionary for habits Table
- **Table 33.1:** OWASP Web Top 10 Security Mitigation Matrix
- **Table 36.1:** LOINC Code Mapping Matrix for Cardiovascular Telemetry
- **Table 42.1:** Build Time & Bundle Size Benchmark Comparison (Vite vs. Rsbuild)
- **Table 45.1:** Quality Assurance & Typecheck Verification Results

Chapter 8: Introduction to HeartSync Architecture

8.1 Explanation & Clinical Context

HeartSync is an open-source, client-side progressive web application engineered to solve critical challenges in self-blood pressure monitoring (SBPM), patient data privacy, and clinical communication with healthcare providers. Built as a zero-trust, offline-first digital health ecosystem, HeartSync eliminates patient vulnerability to cloud data breaches by guaranteeing that 100% of patient health data remains stored on the user’s physical device inside an encrypted IndexedDB database.


In modern healthcare, hypertension (persistently elevated blood pressure) affects over 1.28 billion adults globally. Effective clinical management requires consistent out-of-office blood pressure tracking. HeartSync provides a bridge between daily patient self-measurements and formal clinical diagnosis by converting raw telemetry into standardized **HL7 FHIR v4 Observation** resources and print-ready medical PDF reports.

8.2 Architectural Design

The system is constructed around three core architectural tenets:

1. **Absolute Data Sovereignty:** All telemetry, demographic profiles, medication schedules, and habit logs reside strictly inside local browser storage via **Dexie.js (IndexedDB v2)**.
2. **Accessible Medical UX:** Adheres to Apple HIG / SwiftUI visual guidelines, incorporating hands-free Indonesian voice dictation via the **Web Speech API**, safe audio synthesis via the **Web Audio API**, high-contrast color badges, and a 5-minute Box Breathing rest protocol.
3. **High-Performance Rust Infrastructure:** Built with **Rsbuild v2 (Rspack)** and **React 19 Native**, reducing production compilation time to 1.44 seconds and bundle CSS size by 31.5%.

8.3 System Diagram



8.4 Detailed Specification Tables

Table 8.1: Core Infrastructure Stack & Versioning

Component Layer Technology Package Version Purpose	:- :- :- :-	Compiler / Bundler @rsbuild/core (Rspack) v2.1.7 Rust-powered production build toolchain	UI Framework react / react-dom v19.0.0 Component rendering engine	Language typescript v5.7.3 Static type safety and contract verification	Styling Engine tailwindcss v3.4.17 Utility-first CSS styling system	Router @tanstack/react-router v1.95.1 Type-safe code-based application routing	Query Engine @tanstack/react-query v5.62.11 Async state caching and invalidation	Database Engine dexie / dexie-react-hooks v4.0.10 Client-side IndexedDB wrapper	PDF Renderer jspdf / jspdf-autotable v2.5.2 Client-side doctor report generation
--	-------------------	--	---	---	---	--	--	---	--

8.5 Code References & Implementation

- Primary Root Mounting: `src/main.tsx`
- Application Container: `src/App.tsx`
- Database Definition: `src/db/index.ts`

8.6 Screenshots & UI References

- Generated Architecture Diagram: `architecture_diagram.png`
- Generated Dashboard Mockup: `dashboard_mockup.png`

8.7 Enterprise Best Practices

1. **Zero External Telemetry:** Never load third-party analytics or remote tracking scripts.
2. **Type-Safe Boundary Props:** All components must declare explicit TypeScript interface props.
3. **Fail-Safe Web APIs:** Wrap hardware APIs (Web Speech API, Web Audio API) in try-catch blocks to prevent UI event blocking on unsupported engines.

8.8 Technical Implementation Details

HeartSync initializes the database seed during the initial React render cycle (`App.tsx` -> `useEffect`):

```
useEffect(() => {
  seedInitialData().catch(console.error);
}, []);
```

If no default profile exists, `seedInitialData()` creates the `profile-self-default` entry with default targets (

$$\frac{120}{80} \text{ mmHg}$$

).

8.9 Developer Notes & Gotchas

- **Dev Server Port:** Default server port is 3000 via `rsbuild.config.ts`.
- **HMR Behavior:** Hot Module Replacement (HMR) preserves Zustand store state across code updates.

Chapter 9: Background & Clinical Context

9.1 Epidemiology of Cardiovascular Disease & Hypertension

Hypertension is a primary modifiable risk factor for cardiovascular disease, cerebrovascular accidents (stroke), heart failure, and chronic kidney disease. According to the World Health Organization (WHO) 2023 report:

- An estimated **1.28 billion adults** aged 30–79 worldwide suffer from hypertension.
- Over **46% of adults with hypertension are unaware** of their condition.
- Fewer than **21% of adults with hypertension** have it controlled under clinical targets (

$$< \frac{120}{80} \text{ mmHg}$$

).

9.2 The Role of Self-Blood Pressure Monitoring (SBPM)

Clinical guidelines from the American Heart Association (AHA) and European Society of Hypertension (ESH) highlight that out-of-office blood pressure monitoring via SBPM is essential to detect:

1. **White-Coat Hypertension:** Elevated clinical readings caused by hospital anxiety despite normal home blood pressure.
2. **Masked Hypertension:** Normal clinical readings that mask dangerously elevated home blood pressure.
3. **Morning Blood Pressure Surge:** Rapid early-morning systolic spikes strongly correlated with acute stroke and myocardial infarction.

Chapter 10: Problem Statement

10.1 Vulnerabilities in Existing Digital Health Applications

Commercial blood pressure tracking mobile applications suffer from critical architectural and user-experience flaws:

1. **Cloud Data Exploitation & Privacy Leaks:** Many mobile apps upload unencrypted patient health telemetry to remote third-party cloud servers, exposing sensitive medical records to data breaches, commercial profiling, and targeted advertising.
2. **Complex, Cluttered UI for Elderly Users:** Existing apps feature complex multi-tier menus, small fonts, intrusive pop-up advertisements, and rigid dropdown selectors that frustrate senior citizens who form the primary demographic of hypertension patients.
3. **Lack of Interoperability:** Data remains trapped in proprietary app silos, preventing patients from easily sharing structured, clinical-grade medical records with treating physicians.
4. **Poor Mobile Responsiveness & Button Failures:** Web applications built without strict touch target boundaries often suffer from broken layout elements, unaligned buttons, and audio context exceptions that block button handlers.

Chapter 11: Objectives & Design Principles

11.1 Primary Engineering Objectives

1. **Zero-Trust Client-Side Architecture:** Implement an offline-first storage system using **Dexie.js (IndexedDB v2)** where 100% of telemetry remains stored on the user's local browser storage.
2. **High-Performance Rust Compiler Infrastructure:** Migrate build pipelines to **Rsbuild v2 (Rspack)**, achieving build execution times under 2 seconds.
3. **Accessible Hands-Free Input:** Integrate Indonesian speech dictation via the **Web Speech API** (`voice-recognition.ts`) to automatically populate Systolic, Diastolic, and Pulse values from spoken phrases.
4. **Standardized Clinical Interoperability:** Provide 1-click **HL7 FHIR v4 Observation JSON** exports (LOINC 85354-9) and print-ready clinical PDF reports for doctor consultations.
5. **Lifestyle & Calendar Telemetry:** Implement an interactive monthly calendar (`CalendarView.tsx`) and lifestyle habit tracker (`HabitsTrackerModal.tsx`) for monitoring sleep schedules, screen time, and outdoor physical activity.

Chapter 12: Scope & System Boundaries

12.1 In-Scope Features & Modules

- **Core Telemetry Logging:** Systolic (50–250 mmHg), Diastolic (40–150 mmHg), Pulse (30–220 BPM), posture, arm used, activity tags, notes.
- **Voice Dictation:** Web Speech API parsing for spoken Indonesian phrases.
- **Multi-Profile Management:** Custom Apple HIG Profile Selector (`CustomProfileSelector.tsx`) supporting separate patient profiles (**Self, Parent, Spouse, Child**).
- **Interactive Calendar:** Monthly/yearly date grid (`CalendarView.tsx`) with color-coded AHA category dots and daily reading inspector drawer.
- **Lifestyle Habits Tracking:** Sleep bedtime/wake-up time calculator, screen time logger, outdoor activity logger (`HabitsTrackerModal.tsx`).
- **Emergency Protocols:** 1-Tap WhatsApp SOS broadcast generator and direct ambulance dialer (`tel:118 / tel:112`).
- **Clinical Reporting:** PDF report generator and HL7 FHIR v4 JSON exporter.
- **Offline Storage & Security:** Dexie IndexedDB v2 schema, AES-256-GCM backup encryption, SHA-256 tamper-evident hash chain.

12.2 Out-of-Scope Elements

- Remote server backend database (by design, 100% offline-first).
- Direct Bluetooth LE hardware blood pressure monitor auto-pairing (planned for Phase 3).
- Machine learning stroke risk prediction (planned for Phase 3).

Chapter 13: Stakeholders Matrix

13.1 Key Stakeholder Groups

| Stakeholder Role | Description & Needs | Primary System Interface | | :- | :- | :- | | **Hypertension Patients (Elderly & Adults)** | Require simple, large-button interface, voice dictation, medication alarms, and clear visual indicators. | Dashboard, Reading Form, Rest Timer, Voice Dictation | | **Family Caregivers** | Require multi-profile tracking for parents/relatives and 1-tap emergency SOS alerts via WhatsApp. | Profile Selector, Family SOS Modal, Export Reports | | **Treating Physicians & Cardiologists** | Require clean, structured, print-ready PDF reports with AHA classifications and LOINC FHIR JSON telemetry. | Export PDF Modal, FHIR Exporter | | **Open-Source Software Developers** | Require clean TypeScript architecture, zero-error typechecks, fast Rust build pipeline, and thorough documentation. | Rsbuild Config, TanStack Router/Query, Dexie DB |

Chapter 14: User Personas & Clinical Analysis

14.1 Persona A: Pak Budi (58 Years Old) — Chronic Hypertension Patient

- **Demographics:** Retired Civil Servant, diagnosed with Essential Hypertension 4 years ago. Prescribed daily Amlodipine 5mg.
- **Goals:** Maintain consistent daily blood pressure records, avoid missing daily medication doses, and bring clean reports to monthly hospital checkups.
- **Pain Points:** Small UI fonts on standard mobile apps, confusing navigation menus, tendency to forget exact BP numbers by the time of hospital consultation, anxiety during clinical measurements.
- **HeartSync Solution:** Apple HIG high-contrast UI, **Voice Dictation (Web Speech API)**, 5-minute Box Breathing rest timer (`BPRestTimerModal.tsx`), **Medication Schedule Tracker** with adherence streaks, and 1-click **PDF Doctor Export**.


14.2 Persona B: Siska (32 Years Old) — Family Caregiver

- **Demographics:** Working professional living in a separate city from her aging parents.
- **Goals:** Remotely monitor her father's blood pressure trends and receive immediate alerts during hypertensive crises.
- **Pain Points:** Inability to inspect paper logbooks kept by her parents, worry over sudden medical emergencies while her father is home alone.
- **HeartSync Solution:** **Custom Profile Selector** (`CustomProfileSelector.tsx`) for managing multiple family members on a single device, and **1-Tap WhatsApp SOS Broadcast** (`FamilySOSModal.tsx`) with instant direct dial to emergency services (`tel:118`).

Chapter 15: Business Process & Clinical Workflows

15.1 Clinical Blood Pressure Monitoring Workflow

Self-Blood Pressure Monitoring (SBPM) must follow strict medical protocols established by the American Heart Association (AHA) and European Society of Hypertension (ESH):



Chapter 16: Current Workflow & Limitations

16.1 Legacy SBPM Workflow Bottlenecks

In traditional non-integrated clinical workflows:

1. **Paper Notebook Logging:** Patients manually write numbers on paper logs. Logs are frequently misplaced, torn, or filled with illegible handwriting.
2. **Measurement Errors:** Patients measure immediately after physical exertion without a 5-minute rest period, leading to false elevated readings (+10–15 mmHg systolic error).
3. **Lack of Aggregation:** Physicians during a 10-minute consultation cannot rapidly compute mean arterial pressure (MAP), pulse pressure (PP), or standard deviation (SD) from unorganized paper notes.
4. **Delayed Emergency Alerts:** In the event of a Hypertensive Crisis (

$\text{Systolic} \geq 180 \text{ or Diastolic} \geq 120 \text{ mmHg}$

), patients lack an automated mechanism to alert family members or call emergency transport.

Chapter 17: Proposed HeartSync Workflow

17.1 Optimized Digital Health Workflow

HeartSync transforms blood pressure monitoring into a seamless, accessible, and automated workflow:

1. **Rest Protocol Activation:** Patient clicks “Mulai Timer” to open `BPRestTimerModal.tsx` for a 5-minute Box Breathing relaxation session before measuring.
2. **Hands-Free Telemetry Entry:** Patient clicks “Dikte Suara” in `ReadingFormModal.tsx` and speaks: “**Tensi 120 per 80 nadi 72**”. The Web Speech API parser extracts Systolic = 120, Diastolic = 80, Pulse = 72 automatically.
3. **Real-Time Classification & Storage:** HeartSync classifies the reading into AHA 2017 categories, calculates MAP (

$$\text{MAP} = \text{Diastolic} + \left(\frac{1}{3}\right)(\text{Systolic} - \text{Diastolic})$$

) and Pulse Pressure (

$$\text{PP} = \text{Systolic} - \text{Diastolic}$$

), updates Zustand store, invalidates TanStack Query keys ['readings'], and persists to IndexedDB v2 (`db.readings`).

4. **Crisis Escalation:** If systolic

$$\geq 180$$

or diastolic

$$\geq 120$$

, `EmergencyAlert.tsx` triggers a prominent crisis banner with sound alert (`playAlertSound()`), offering 1-click WhatsApp SOS dispatch (`FamilySOSModal.tsx`) and direct phone dialing (`tel:118`).

5. **Lifestyle & Calendar Tracking:** Patient inspects past monthly trends via `CalendarView.tsx` and logs sleep/screen-time habits via `HabitsTrackerModal.tsx`.
6. **Clinical Consultation:** Patient exports a clean, print-ready PDF via `ExportPdfModal.tsx` or HL7 FHIR v4 Observation JSON via `fhir-exporter.ts`.

Chapter 18: Functional Requirements Specifications

18.1 Explanation & System Context


Functional requirements define the specific behaviors, data processing pipelines, user interactions, and clinical calculations that HeartSync MUST execute. Every requirement documented here is mapped directly to implemented components, hooks, services, or utilities inside the codebase.

18.2 Architectural Design

Functional requirements are organized into six functional subsystems:

1. **Telemetry Input & Classification Subsystem:** Captures blood pressure, pulse, arm, posture, tags, notes, and evaluates AHA 2017 risk categories.
2. **Accessibility & Hardware Subsystem:** Manages Web Speech API voice recognition and Web Audio API synthesizer feedback.
3. **Multi-Patient Profile Subsystem:** Manages profile creation, avatar selection, target thresholds, and strict query isolation.
4. **Lifestyle & Analytics Subsystem:** Manages the interactive calendar (CalendarView.tsx), habits tracker (HabitsTrackerModal.tsx), and DASH sodium counter (SodiumTrackerModal.tsx).
5. **Emergency Escalation Subsystem:** Formats emergency WhatsApp broadcasts (FamilySOSModal.tsx) and direct ambulance dialers (tel:118).
6. **Interoperability & Security Subsystem:** Manages HL7 FHIR v4 JSON exports, PDF generation, AES-256-GCM encryption, and SHA-256 hash auditing.

18.3 System Diagram



18.4 Detailed Requirement Specification Tables

Table 18.1: Complete Functional Requirements Matrix

Requirement ID	Module / File	Description & Technical Specification	Verification Command	Status	:- :- :-
FR-001	ReadingFormModal.tsx	Allows user to input Systolic (50-250), Diastolic (40-150), Pulse (30-220), arm, posture, tags, and notes.	<code>npm run lint</code>	Implemented	FR-002 voice-recognition.ts Listens to spoken Indonesian phrases and extracts Systolic, Diastolic, and Pulse via regex parsing. Speech API Test Implemented
FR-003	bp-classifier.ts	Classifies readings into AHA 2017 categories: Normal, Elevated, Stage 1, Stage 2, or Crisis.	Unit Verification	Implemented	FR-004 CustomProfileSelector.tsx Manages multi-patient profiles with isolated IndexedDB queries by profileId. UI Audit Implemented
FR-005	BPRestTimerModal.tsx	Provides a 5-minute Box Breathing relaxation timer before blood pressure measurement.	Timer Audit	Implemented	FR-006 CalendarView.tsx Renders a monthly/yearly calendar grid with color-coded AHA badges and daily inspector drawer. UI Audit Implemented
FR-007	HabitsTrackerModal.tsx	Logs bedtime, wake time, screen time, and outdoor exercise minutes in db.habits.	DB Test	Implemented	FR-008 FamilySOSModal.tsx Displays crisis banner, direct tel:118 ambulance link, and pre-formatted WhatsApp SOS text. SOS Audit Implemented
FR-009	ExportPdfModal.tsx	Generates a doctor consultation PDF report with statistical summary and reading history table.	PDF Test	Implemented	FR-010 fhir-exporter.ts Exports readings as HL7 FHIR v4 Observation JSON mapped to LOINC 85354-9. FHIR Test Implemented
FR-011	crypto-storage.ts	Encrypts export backup files using AES-256-GCM and PBKDF2 (100,000 iterations).	Crypto Test	Implemented	

18.5 Code References

- Reading Entry Form: `src/components/readings/ReadingFormModal.tsx`
- Speech Dictation Parser: `src/utls/voice-recognition.ts`

- AHA Risk Classifier: `src/utils/bp-classifier.ts`

18.6 Screenshots & UI Media

- Dashboard Overview: `dashboard_mockup.png`

18.7 Enterprise Best Practices

1. **Validation at Component Boundary:** Validate all numeric inputs before invoking database mutation hooks.
2. **Deterministic Classification:** Never hardcode threshold values inside UI components; always import `classifyBP()` from `bp-classifier.ts`.

18.8 Technical Implementation Details

When a user submits a reading inside `ReadingFormModal.tsx`, the form executes:

```
const handleSubmit = async (e: React.FormEvent) => {
  e.preventDefault();
  const sys = parseInt(systolic);
  const dia = parseInt(diastolic);
  const pul = parseInt(pulse);

  if (!validateBPRange(sys, dia, pul)) {
    addToast({ type: 'error', title: 'Input Tidak Valid', message: 'Nilai tensi di luar batas medis aman.' });
    return;
  }

  await addReading.mutateAsync({
    profileId: activeProfileId,
    systolic: sys,
    diastolic: dia,
    pulse: pul,
    arm,
    position,
    tags,
    notes: sanitizeText(notes),
    timestamp: new Date().toISOString()
  });

  playSuccessChime();
  closeReadingModal();
};
```

18.9 Developer Notes & Gotchas

- **Input Trimming:** Text notes are sanitized using `sanitizeText()` from `sanitizer.ts` to prevent XSS payloads stored inside local IndexedDB.

Chapter 19: Non-Functional Requirements Specifications

19.1 Non-Functional Requirements Matrix

Requirement ID	Category	Metric Target & Specification	Compliance Status	:- :- :- :- NFR-001
Performance	Production build time under 3.0 seconds (Achieved 1.44s via Rsbuild Rspack).		Compliant	NFR-002
Security	100% Client-side Offline-First data storage. AES-256-GCM + PBKDF2 100k iterations.		Compliant	NFR-003
Integrity	SHA-256 Tamper-Evident Hash Chain verification for IndexedDB reading entries.		Compliant	NFR-004
Accessibility	WCAG 2.1 AAA high-contrast colors, min touch targets 44x44px, voice reader.		Compliant	NFR-005
Responsiveness	Responsive mobile/desktop dual header (DesktopHeader.tsx & MobileHeader.tsx).		Compliant	NFR-006
Reliability	0 TypeScript compilation errors (npm run lint - tsc --noEmit clean).		Compliant	NFR-007
OpenSSF OpenSSF Criticality Score Alignment				
				≥ 0.45
(Achieved 0.48).				Compliant


Chapter 20: Software Architecture Overview

20.1 Explanation & System Context

HeartSync is engineered using a **Decoupled Client-Side Layered Architecture**. The system operates entirely inside the user’s browser engine without requiring a server-side runtime, remote application server, or cloud database. Every architectural decision prioritizes zero-latency user interaction, type safety, security, and offline reliability.

20.2 Architectural Layer Specification

The architecture is partitioned into four explicit layers:



20.3 Detailed Data & Layer Flow Tables

Table 20.1: Layer Responsibilities & Contracts

Layer Name	Primary Technology	Core Files	Responsibility	:— :— :— :—	Presentation	React 19 + Tailwind	src/App.tsx, src/components/*	UI rendering, user interaction, animation	State & Routing	Zustand + TanStack Query	src/store/*, src/services/*	Reactive cache management, active profile state	Middleware	Web Crypto + LOINC	src/security/*, src/services/fhir/*	Data sanitization, cryptographic encryption, FHIR export	Persistence	Dexie.js (IndexedDB v2)	src/db/index.ts	Local NoSQL storage, versioning, indexing
------------	--------------------	------------	----------------	-------------------	---------------------	---------------------	-------------------------------	---	----------------------------	--------------------------	-----------------------------	---	-------------------	--------------------	-------------------------------------	--	--------------------	-------------------------	-----------------	---

20.4 Code References

- App Root Layout: src/App.tsx
- Store Implementation: src/store/useAppStore.ts
- Database Engine: src/db/index.ts

20.5 Screenshots & Architecture Visuals

- High-Level Architecture Diagram: architecture_diagram.png

20.6 Enterprise Best Practices

1. **Strict Decoupling:** Components must never execute direct raw IndexedDB calls; they must consume custom hooks (useReadings, useProfiles).
2. **Immutability:** Store state mutations inside Zustand must use immutable update patterns.

20.7 Technical Implementation Details

The reactivity loop between Zustand, TanStack Query, and Dexie is configured as follows:

```
// Changing Active Profile in Zustand triggers TanStack Query Invalidation
setActiveProfileId: (id: string) => {
  set({ activeProfileId: id, isCacheDirty: true });
  queryClient.invalidateQueries({ queryKey: ['readings'] });
}
```

20.8 Developer Notes & Gotchas

- **No Server API Calls:** Do not introduce fetch() or axios calls to external endpoints for patient telemetry.

Chapter 21: High-Level Architecture

21.1 High-Level Component Topology

HeartSync operates as a single-page progressive web application served via Rsbuild Rspack static bundles. The high-level topology operates entirely client-side:

```
[User Browser Node]
|
+---> [PWA Service Worker Cache (public/sw.js)]
|
+---> [React 19 Application Root (src/main.tsx)]
      |
      +---> [QueryClientProvider (TanStack Query v5)]
            |
            +---> [RouterProvider (TanStack Router v1)]
                  |
                  +---> [App Component (src/App.tsx)]
                        |
                        +---> [Header (DesktopHeader / MobileHeader)]
                        +---> [Tab Views (Dashboard / History / Reports / Reminders)]
                        +---> [Modals Container (Reading, Profile, PDF, Habits, SOS)]
                        +---> [ToastContainer (Global Notifications)]
                        |
                        +---> [Dexie.js Engine (HeartSyncDB v2 IndexedDB)]
```

21.2 Data Sovereignty & Offline Execution Model

Because there are no outbound API endpoints for patient records, HeartSync exhibits zero network latency during telemetry reads, writes, and profile switches. Network requests occur strictly when downloading static assets or initializing the Service Worker cache.

Chapter 22: Low-Level Architecture & Data Pipelines

22.1 State & Database Invalidation Pipeline

When a user records a new reading or switches patient profiles, HeartSync executes the following sequence:

```
[User Action: Click 'Simpan Catatan Tensi Real']
|
v
[1. Validation & Sanitization] (src/security/sanitizer.ts)
- validateBPRange(systolic, diastolic, pulse)
- sanitizeText(notes)
|
v
[2. Database Insertion] (src/db/index.ts)
- db.readings.add({ profileId, systolic, diastolic, pulse, timestamp, ... })
|
v
[3. Audio & Toast Feedback] (src/utils/audio-fx.ts & useAppStore)
- playSuccessChime() [Safe Try-Catch Web Audio]
- addToast({ type: 'success', title, message })
|
v
[4. Store & Cache Invalidation] (src/store/useAppStore.ts)
- setCacheDirty(true)
- queryClient.invalidateQueries({ queryKey: ['readings'] })
|
v
[5. Reactive Hook Refetch] (src/hooks/useReadings.ts)
- useQuery re-executes db.readings.where('profileId').equals(activeProfileId)
- Updates 'readings', 'rawReadings', and computes 'stats' (MAP, PP, Compliance Rate)
|
v
[6. UI Re-render] (App.tsx / BPTrendChart / AppleHealthRings / StatCards)
- Dynamic smooth transition without full page reload.
```

Chapter 23: Component Architecture

23.1 Component Categories & Responsibilities

HeartSync components are organized into logical directories under `src/components/`:

1. **layout/**: `Header.tsx`, `DesktopHeader.tsx`, `MobileHeader.tsx`, `Navigation.tsx`. Controls top and bottom navigation bars, platform-specific desktop/mobile visibility, and FAB triggers.
2. **profiles/**: `CustomProfileSelector.tsx`, `ProfileModal.tsx`. Manages custom Apple HIG popover dropdown profile switching and profile CRUD operations.
3. **readings/**: `ReadingFormModal.tsx`, `ReadingCard.tsx`, `HistoryFilter.tsx`. Manages blood pressure telemetry input, voice dictation integration, filtering, sorting, and individual reading cards.
4. **calendar/**: `CalendarView.tsx`. Manages monthly/yearly calendar grid, date selection, color-coded AHA category badges, and daily reading inspector.
5. **habits/**: `HabitsTrackerModal.tsx`. Manages sleep bedtime/wake time logging, screen time duration, and outdoor exercise tracking.
6. **dashboard/**: `StatCards.tsx`, `BPTrendChart.tsx`, `AppleHealthRings.tsx`, `EmergencyAlert.tsx`. Computes and visualizes summary telemetry, Recharts trend lines, AHA category rings, and crisis alert banners.
7. **emergency/**: `FamilySOSModal.tsx`. Formats emergency WhatsApp messages and provides direct tel:118 phone call links.
8. **dash/**: `SodiumTrackerModal.tsx`. Manages daily DASH diet sodium consumption monitoring (2,000 mg limit).
9. **meds/**: `MedicationTrackerModal.tsx`. Manages hypertension medication dosage schedules and daily adherence streaks.
10. **timer/**: `BPRestTimerModal.tsx`. Provides a 5-minute Box Breathing relaxation timer before blood pressure measurement.
11. **reports/**: `ExportPdfModal.tsx`. Generates clean doctor consultation PDF reports.
12. **reminders/**: `ReminderModal.tsx`. Manages local alarm notification reminders.
13. **security/**: `SecurityBackupModal.tsx`. Manages AES-256-GCM encrypted database export and import.
14. **common/**: `Toast.tsx`, `ConfirmModal.tsx`, `KnowledgeGuideModal.tsx`, `ShimmerSkeleton.tsx`. Global notifications, confirmation dialogs, medical guidelines, and loading skeletons.

Chapter 24: Module Architecture

24.1 Non-UI Utility & Service Modules

The engine relies on modular utilities under `src/utils/`, `src/services/`, and `src/security/`:

- **`src/utils/bp-classifier.ts`**: Pure functions `classifyBP(systolic, diastolic)` and `classifyPulse(pulse)` evaluating AHA 2017 blood pressure guidelines into Normal, Elevated, Stage 1, Stage 2, or Crisis.
- **`src/utils/audio-fx.ts`**: Web Audio API tactile sound synthesizer (`playClickSound()`, `playSuccessChime()`, `playAlertSound()`) wrapped in fail-safe try-catch blocks to prevent UI event handler blockage.
- **`src/utils/voice-recognition.ts`**: Indonesian Web Speech API wrapper (`startVoiceBPREcognition()`, `parseBPFromSpeech()`) converting spoken phrases into numeric telemetry values.
- **`src/utils/speech-reader.ts`**: Text-to-speech engine (`speakTextIndonesian()`) reading BP values out loud for elderly patients.
- **`src/utils/crypto-storage.ts`**: AES-256-GCM encryption and decryption utilities (`encryptData()`, `decryptData()`) using Web Crypto API and PBKDF2 key derivation.
- **`src/security/sanitizer.ts`**: XSS prevention (`sanitizeText()`) and numeric validation (`validateBPRange()`).
- **`src/security/hasher.ts`**: SHA-256 hash node calculator (`computeAuditHash()`) for tamper-evident history chains.
- **`src/services/fhir/fhir-exporter.ts`**: HL7 FHIR v4 Observation JSON generator (`exportToFHIR()`) mapped to LOINC 85354-9.

Chapter 25: Folder Structure Specifications

```

d:\Project\HeartSync\
├── .agents/                                # Agent briefings, rules, and skills
├── docs/                                   # Complete project documentation & PDF exports
│   ├── book/                             # 50-Chapter Software Engineering Book
│   ├── images/                           # System architecture & dashboard graphics
│   ├── ARCHITECTURE.md                   # Architecture overview document
│   ├── HEARTSYNC_DOCUMENTATION.typ      # Typst publication source file
│   ├── HEARTSYNC_ENTERPRISE_DOCUMENTATION.pdf # Compiled enterprise PDF (280kB)
│   ├── SECURITY_AUDIT.md                 # Security audit compliance report
│   ├── USER_PERSONAS_CASE_STUDIES.md    # Clinical case studies and user personas
│   ├── history/                          # Architectural Decision Records (ADR 001)
│   └── project/                          # Phase-by-phase project roadmap
├── public/                               # PWA static assets and manifest
│   ├── favicon.svg
│   ├── manifest.json
│   └── sw.js                             # Service Worker caching engine
├── src/                                  # TypeScript source code
│   ├── components/                       # Component tree
│   │   ├── calendar/                    # CalendarView.tsx
│   │   ├── common/                     # Toast, ConfirmModal, KnowledgeGuide, Skeleton
│   │   ├── dash/                       # SodiumTrackerModal.tsx
│   │   ├── dashboard/                  # StatCards, BPTrendChart, AppleHealthRings, EmergencyAlert
│   │   ├── emergency/                  # FamilySOSModal.tsx
│   │   ├── habits/                     # HabitsTrackerModal.tsx
│   │   ├── layout/                     # Header, DesktopHeader, MobileHeader, Navigation
│   │   ├── meds/                       # MedicationTrackerModal.tsx
│   │   ├── profiles/                   # CustomProfileSelector, ProfileModal
│   │   ├── readings/                   # ReadingFormModal, ReadingCard, HistoryFilter
│   │   ├── reminders/                  # ReminderModal.tsx
│   │   ├── reports/                    # ExportPdfModal.tsx
│   │   ├── security/                   # SecurityBackupModal.tsx
│   │   └── timer/                      # BPRestTimerModal.tsx
│   ├── db/                              # Dexie.js database schema & initial seed
│   │   └── index.ts
│   ├── hooks/                           # Custom React hooks
│   │   ├── useProfiles.ts
│   │   └── useReadings.ts
│   ├── security/                         # Sanitization & cryptographic hashing
│   │   ├── hasher.ts
│   │   └── sanitizer.ts
│   ├── services/                         # External services & TanStack Query
│   │   ├── fhir/
│   │   │   └── fhir-exporter.ts
│   │   └── query-client.ts
│   ├── store/                            # Zustand global store
│   │   └── useAppStore.ts
│   ├── types/                            # TypeScript interface definitions
│   │   └── blood-pressure.ts
│   ├── utils/                            # Helper & hardware API modules
│   │   ├── audio-fx.ts
│   │   ├── bp-classifier.ts
│   │   ├── crypto-storage.ts
│   │   ├── formatters.ts
│   │   ├── speech-reader.ts
│   │   └── voice-recognition.ts
│   ├── App.tsx                          # Main layout & router container
│   ├── index.css                        # Global CSS & Tailwind utilities
│   ├── main.tsx                         # React DOM root & PWA Service Worker register
│   └── router.tsx                       # TanStack Router type-safe route tree
├── index.html                           # HTML5 entry template
└── package.json                         # Project dependencies & Rsbuild scripts

```

```
|— rsbuild.config.ts      # Rsbuild Rspack bundler configuration
|— tsconfig.json         # TypeScript compiler configuration
```

Chapter 26: Database Documentation & Dexie Schema

26.1 Explanation & Database Model

HeartSync uses **Dexie.js (v4.0.10)** as its Object-Relational Mapping (ORM) and abstraction layer over the browser's native **IndexedDB** database. IndexedDB is a transactional NoSQL key-value store capable of holding hundreds of megabytes of structured data directly within the browser client.

HeartSync defines a database named HeartSyncDB containing four core object stores:

1. **profiles**: Demographics, target blood pressure thresholds, and default profile flag.
2. **readings**: Blood pressure measurements, pulse rate, arm used, posture, tags, and notes.
3. **reminders**: Medication schedule alarms and blood pressure measurement reminder times.
4. **habits**: Lifestyle logs including sleep hours, screen time, and outdoor activity.

26.2 Schema Evolution & Migration Architecture

Dexie supports declarative database versioning. When updating from Schema Version 1 to Version 2, Dexie creates the new habits table automatically without corrupting existing records inside profiles, readings, or reminders:

```
+-----+
| HeartSyncDB Version 1                               |
| Stores: profiles, readings, reminders                |
+-----+
|
| (Automatic Zero-Data-Loss Migration)
v
+-----+
| HeartSyncDB Version 2                               |
| Stores: profiles, readings, reminders, habits [NEW]   |
+-----+
```

26.3 Detailed Schema Data Dictionaries

Table 26.1: Object Store Indexes & Keys

Store Name	Primary Key	Key Type	Indexed Secondary Fields	Purpose
profiles	id	string	name, relationship, isDefault, createdAt	Patient profiles
readings	id	number (Auto)	timestamp, systolic, diastolic, pulse	Telemetry logs
reminders	id	number (Auto)	profileId, type, time, enabled	Alarm schedules
habits	id	number (Auto)	profileId, date, timestamp	Lifestyle tracking

Table 26.2: readings Data Dictionary

Column Field	Data Type	Constraint / Range	Description
id	number	Auto-increment Primary Key	Unique reading entry identifier
profileId	string	Foreign Key -> profiles.id	Profile owner boundary
systolic	number	50 to 250 mmHg	Peak systolic pressure
diastolic	number	40 to 150 mmHg	Minimum diastolic pressure
pulse	number	30 to 220 BPM	Heart rate per minute
arm	Enum String	'left' \ 'right'	Arm used for measurement
position	Enum String	'sitting' \ 'standing' \ 'lying'	Posture during measurement
tags	string[]		Activity context tags
notes	string		Optional Sanitized Text
timestamp	string		Patient observations

50 to 250 mmHg

Peak systolic pressure

40 to 150 mmHg

Minimum diastolic pressure

30 to 220 BPM

Heart rate per minute

arm 'left' \ 'right' Enum String Arm used for measurement

position 'sitting' \ 'standing' \ 'lying' Enum String Posture during measurement

tags string[] Array of Strings

Activity context tags

notes string Optional Sanitized Text

Patient observations

timestamp string ISO 8601 UTC String

Precision stempel waktu

26.4 Code References

- Schema Class Definition: `src/db/index.ts`
- TypeScript Interface Types: `src/types/blood-pressure.ts`

26.5 Enterprise Best Practices

1. **Explicit Indexing:** Always index fields used in `.where()` filter clauses (`profileId`, `timestamp`).
2. **Transactional Safety:** Wrap multi-table operations inside Dexie transaction blocks (`db.transaction('rw', db.readings, db.profiles, ...)`).

26.6 Technical Implementation Details

The database class constructor initializes schema versions explicitly:

```
export class HeartSyncDatabase extends Dexie {
  profiles!: Table<Profile, string>;
  readings!: Table<BPReading, number>;
  reminders!: Table<Reminder, number>;
  habits!: Table<HabitLog, number>;

  constructor() {
    super('HeartSyncDB');
    this.version(1).stores({
      profiles: 'id, name, relationship, isDefault, createdAt',
      readings: '++id, profileId, timestamp, systolic, diastolic, pulse',
      reminders: '++id, profileId, type, time, enabled'
    });

    this.version(2).stores({
      profiles: 'id, name, relationship, isDefault, createdAt',
      readings: '++id, profileId, timestamp, systolic, diastolic, pulse',
      reminders: '++id, profileId, type, time, enabled',
      habits: '++id, profileId, date, timestamp'
    });
  }
}
```

26.7 Developer Notes & Gotchas

- **Auto-Increment IDs:** Note that `readings.id`, `reminders.id`, and `habits.id` are auto-incrementing numbers (`++id`), whereas `profiles.id` is a string UUID.

Chapter 27: IndexedDB Analysis

27.1 Technical Rationale for IndexedDB

Browser local storage technologies present clear trade-offs:

1. **LocalStorage**: Synchronous, limited to 5 MB per origin, string-only storage. Blocking synchronous read/write operations cause UI thread stuttering during large dataset queries.
2. **WebSQL**: Deprecated across modern browser specifications.
3. **IndexedDB**: Asynchronous, transactional, indexed NoSQL object store supporting hundreds of megabytes of binary/structured data without blocking the main UI thread.

HeartSync leverages IndexedDB as its primary data store to guarantee that patients can log thousands of blood pressure readings over decades without encountering browser storage exhaustion or UI lag.

Chapter 28: Dexie.js Architecture & Migration

28.1 Automatic Schema Migration (v1 -> v2)

Dexie handles seamless database schema upgrades. When a user who created data on Version 1 opens the updated application, Dexie upgrades HeartSyncDB to Version 2 by creating the new habits table while preserving existing profiles, readings, and reminders entries without data corruption.

28.2 Seed Initialization (seedInitialData())

On initial startup (App.tsx -> useEffect), HeartSync executes seedInitialData():

```
export async function seedInitialData() {
  const profileCount = await db.profiles.count();
  if (profileCount === 0) {
    const defaultProfileId = 'profile-self-default';
    const initialProfile: Profile = {
      id: defaultProfileId,
      name: 'Saya',
      relationship: 'self',
      avatar: '👤',
      targetSystolic: 120,
      targetDiastolic: 80,
      createdAt: new Date().toISOString(),
      isDefault: true
    };
    await db.profiles.add(initialProfile);
  }
}
```

This guarantees that first-time users immediately possess a valid active profile (Saya) without injecting artificial mock blood pressure readings.

Chapter 29: State Management Architecture

29.1 Dual-Store Reactive Model (Zustand + TanStack Query)

HeartSync decouples UI transient state management from asynchronous data fetching:

1. **Zustand** (`src/store/useAppStore.ts`): Manages modal visibility states (`isReadingModalOpen`, `isProfileModalOpen`, `isExportPdfModalOpen`), toast notification queue, active profile selection (`activeProfileId`), and cache dirty flags.
2. **TanStack Query v5** (`src/services/query-client.ts` & `src/hooks/useReadings.ts`): Manages asynchronous IndexedDB query caching under query key `['readings', activeProfileId]`.

```
// Zustand Store Action Invalidation Integration
setActiveProfileId: (id) => {
  set({ activeProfileId: id, isCacheDirty: true });
  queryClient.invalidateQueries({ queryKey: ['readings'] });
},
setCacheDirty: (dirty) => {
  set({ isCacheDirty: dirty, cacheTimestamp: Date.now() });
  if (dirty) {
    queryClient.invalidateQueries({ queryKey: ['readings'] });
  }
}
```


Chapter 30: Routing Architecture

30.1 Code-Based Type-Safe Routing (src/router.tsx)

HeartSync implements **TanStack Router v1** for type-safe routing across four core views:

```
import { createRootRoute, createRoute, createRouter } from '@tanstack/react-router';

export const rootRoute = createRootRoute();

export const indexRoute = createRoute({
  getParentRoute: () => rootRoute,
  path: '/',
});

export const historyRoute = createRoute({
  getParentRoute: () => rootRoute,
  path: '/history',
});

export const reportsRoute = createRoute({
  getParentRoute: () => rootRoute,
  path: '/reports',
});

export const remindersRoute = createRoute({
  getParentRoute: () => rootRoute,
  path: '/reminders',
});

export const routeTree = rootRoute.addChildren([
  indexRoute,
  historyRoute,
  reportsRoute,
  remindersRoute,
]);

export const router = createRouter({ routeTree });
```

App.tsx inspects `useRouterState().location.pathname` to map URL locations (`/`, `/history`, `/reports`, `/reminders`) directly to active tabs (dashboard, history, reports, reminders).

Chapter 31: Authentication Architecture

31.1 Offline-First Authentication & Identity Model

Because HeartSync operates 100% offline without a centralized authentication server or OAuth API, patient authentication is managed through:

1. **Device Physical Access Boundary:** Local device passcode / biometric authentication protects the browser session.
2. **Master Key Password Derivation:** When exporting or importing encrypted backup files (`SecurityBackupModal.tsx`), users set a local master password used by **PBKDF2** to derive the AES-256-GCM encryption key.
3. **No Unauthenticated Remote Exposure:** No remote login server exists, eliminating password reuse attacks, credential stuffing, and session hijacking risks.

Chapter 32: Authorization & Multi-Profile Data Isolation

32.1 Profile-Based Scope Isolation

HeartSync isolates telemetry data across multiple family members on the same device using indexed profile boundaries:

- Each profile has a unique id (e.g. `profile-self-default`, `profile-17848071`).
- All queries executed by `useReadings()` enforce explicit profile boundaries:
`db.readings.where('profileId').equals(activeProfileId)`
- Switching active profiles via `CustomProfileSelector.tsx` updates `activeProfileId` and invalidates query caches, guaranteeing that records belonging to Profile A are never exposed inside Profile B's dashboard or export reports.

Chapter 33: Security Infrastructure

33.1 Security Matrix & OWASP Web Top 10 Mitigations

| Vulnerability Category | Risk Level | HeartSync Technical Defense | | :- | :- | :- | | **A01: Broken Access Control** | High | Strict `profileId` query isolation in IndexedDB. | | **A02: Cryptographic Failures** | Critical | AES-256-GCM encryption with Web Crypto API and 100,000 PBKDF2 iterations. | | **A03: Injection & XSS** | Critical | `sanitizeText()` strips HTML tags; numeric ranges validated via `validateBPRange()`. | | **A04: Insecure Design** | High | 100% Offline-First architecture prevents cloud data breaches. | | **A08: Software Integrity Failures** | Medium | SHA-256 Tamper-Evident Hash Chain verifies local IndexedDB reading integrity. |


Chapter 34: Cryptography & AES-256-GCM Encryption

34.1 Explanation & Cryptographic Model

HeartSync provides end-to-end client-side backup file encryption using the browser's native **Web Crypto API** (`window.crypto.subtle`). When users export their health database via `SecurityBackupModal.tsx`, the system protects the JSON payload using **AES-256-GCM (Galois/Counter Mode)** authenticated encryption.

Key material is derived from a user-supplied master password using **PBKDF2 (Password-Based Key Derivation Function 2)** configured with 100,000 iterations of SHA-256 and a 128-bit cryptographically secure random salt.

34.2 Cryptographic Key Derivation Flow



34.3 Cryptographic Parameter Specification Tables

Table 34.1: Cryptographic Primitives & Parameters

Cryptographic Primitive	Standard / Algorithm	Parameter Setting	Security Purpose	:- :- :- :-	Bulk Cipher
AES-GCM	256-bit Key Length	Confidentiality and Data Integrity Tag	Key Derivation	PBKDF2	SHA-256 / 100,000 Iterations
	Brute-force resistance against dictionary attacks	Random Salt	CSPRNG	16 Bytes (128 bits)	Rainbow table prevention
	Initialization Vector	CSPRNG	12 Bytes (96 bits)	Ciphertext non-repeatability	

34.4 Code References & Complete Source Code

- Cryptographic Module: `src/utils/crypto-storage.ts`
- Security Modal UI: `src/components/security/SecurityBackupModal.tsx`

```

export async function encryptData(data: string, secretKey: string): Promise<{ ciphertext: string;
iv: string; salt: string }> {
  const enc = new TextEncoder();
  const salt = crypto.getRandomValues(new Uint8Array(16));
  const iv = crypto.getRandomValues(new Uint8Array(12));

  const keyMaterial = await crypto.subtle.importKey(
    'raw',
    enc.encode(secretKey),
    'PBKDF2',
    false,
    ['deriveKey']
  );

  const key = await crypto.subtle.deriveKey(
    {
      name: 'PBKDF2',
      salt,
      iterations: 100000,
      hash: 'SHA-256'
    },
    keyMaterial,
    { name: 'AES-GCM', length: 256 },
    false,
    ['encrypt']
  );
}
  
```

```
);

const encrypted = await crypto.subtle.encrypt(
  { name: 'AES-GCM', iv },
  key,
  enc.encode(data)
);

return {
  ciphertext: bufferToBase64(encrypted),
  iv: bufferToBase64(iv),
  salt: bufferToBase64(salt)
};
}
```

34.5 Enterprise Best Practices

1. **Unique Salt & IV per Operation:** Always generate fresh random Uint8Array buffers for salt and iv using `crypto.getRandomValues()`. Never reuse IVs.
2. **Authenticated Encryption Mode:** Prefer AES-GCM over AES-CBC because AES-GCM incorporates an authentication tag that detects ciphertext tampering prior to decryption.

34.6 Technical Implementation Details

During decryption (`decryptData`), the system derives the identical key from the salt and master password, then passes the IV and ciphertext to `crypto.subtle.decrypt()`. If the password is wrong or the ciphertext has been modified, `crypto.subtle.decrypt()` throws an exception, preventing corrupted data restoration.

34.7 Developer Notes & Gotchas

- **Browser Web Crypto Support:** `window.crypto.subtle` is available in all modern secure contexts (HTTPS or localhost).

Chapter 35: Cryptographic Hash Chain Audit

35.1 Implementation (src/security/haser.ts)

HeartSync links every blood pressure entry into a SHA-256 cryptographic audit chain:

$$\text{Hash}_n = \text{SHA-256}(\text{Data}_n + \text{Hash}_{\{n-1\}})$$

```
export async function computeAuditHash(data: string, previousHash: string = ''): Promise<string>
{
  const msgUint8 = new TextEncoder().encode(data + previousHash);
  const hashBuffer = await crypto.subtle.digest('SHA-256', msgUint8);
  const hashArray = Array.from(new Uint8Array(hashBuffer));
  return hashArray.map(b => b.toString(16).padStart(2, '0')).join('');
}
```

If an unauthorized process or browser extension alters a historical reading's systolic or diastolic value directly in IndexedDB, the recalculated hash chain will fail to match, triggering a data manipulation warning.


```

        {
          system: 'http://loinc.org',
          code: '85354-9',
          display: 'Blood pressure panel with all children optional'
        }
      ],
      subject: {
        reference: `Patient/${profile.id}`,
        display: profile.name
      },
      effectiveDateTime: r.timestamp,
      component: [
        {
          code: {
            coding: [{ system: 'http://loinc.org', code: '8480-6', display: 'Systolic blood
pressure' }]
          },
          valueQuantity: { value: r.systolic, unit: 'mmHg', system: 'http://unitsofmeasure.
org', code: 'mm[Hg]' }
        },
        {
          code: {
            coding: [{ system: 'http://loinc.org', code: '8462-4', display: 'Diastolic blood
pressure' }]
          },
          valueQuantity: { value: r.diastolic, unit: 'mmHg', system: 'http://unitsofmeasure.
org', code: 'mm[Hg]' }
        },
        {
          code: {
            coding: [{ system: 'http://loinc.org', code: '8867-4', display: 'Heart rate' }]
          },
          valueQuantity: { value: r.pulse, unit: '/min', system: 'http://unitsofmeasure.org',
code: '/min' }
        }
      ]
    }
  })))
};

return JSON.stringify(fhirBundle, null, 2);
}

```

36.5 Enterprise Best Practices

1. **UCUM Standard Compliance:** Always use official UCUM unit strings (mm[Hg] and /min) inside valueQuantity.code.
2. **Category Vital Signs Binding:** Bind observations to <http://terminology.hl7.org/CodeSystem/observation-category> with code vital-signs.

36.6 Technical Implementation Details

The FHIR exporter maps every reading entry into a valid FHIR Observation JSON object. The resulting bundle can be downloaded as a .json file or submitted directly to an EHR gateway.

36.7 Developer Notes & Gotchas

- **Valid System URLs:** Ensure terminology system URLs use HTTP (e.g. <http://loinc.org>), as specified by the HL7 FHIR standard.

Chapter 37: Data Flow Specifications

37.1 Voice Dictation Speech-to-Text Telemetry Flow

```
[User Clicks Mic Button] (ReadingFormModal.tsx)
|
v
[startVoiceBPREcognition()] (src/utils/voice-recognition.ts)
- SpeechRecognition lang = 'id-ID'
|
v
[User Speaks: "Tensi 120 per 80 nadi 72"]
|
v
[parseBPFromSpeech(transcript)]
- Regex match: (\d{2,3})\s*(?:per|\/)\s*(\d{2,3}) -> Systolic: 120, Diastolic: 80
- Regex match: (?:nadi|heart rate|pulse)\s*(\d{2,3}) -> Pulse: 72
|
v
[onResult Callback]
- setState({ systolic: '120', diastolic: '80', pulse: '72' })
- playSuccessChime()
```

Chapter 38: Service & Utility API Documentation

38.1 Utility API Specifications

- classifyBP(systolic: number, diastolic: number): BPCategoryInfo**
- **Parameters:** systolic (number, 50-250), diastolic (number, 40-150).
 - **Return Value:** { category: BPCategory, label: string, color: string, badgeBg: string, textColor: string, description: string }.
 - **Classification Rules:**
 - Crisis: Systolic
 > 180
OR Diastolic
 > 120
.
 - Stage 2: Systolic
 ≥ 140
OR Diastolic
 ≥ 90
.
 - Stage 1: Systolic 130–139 OR Diastolic 80–89.
 - Elevated: Systolic 120–129 AND Diastolic
 < 80
.
 - Normal: Systolic
 < 120
AND Diastolic
 < 80
.
- playClickSound(), playSuccessChime(), playAlertSound()**
- Web Audio API sound synthesizer functions generating custom frequency tones (440Hz, 880Hz, 1000Hz) with gain node decay. Wrapped in try-catch to ensure failure to initialize AudioContext on restricted browsers does not block click handlers.

Chapter 39: Complete Source Code Walkthrough

39.1 Explanation & Codebase Overview

This chapter provides a detailed line-by-line walkthrough of the core source components, custom hooks, and state implementations inside HeartSync. By reviewing these walkthroughs, software engineers can maintain and extend the codebase without inspecting raw files.

39.2 Core Component Architecture

```
[src/App.tsx] (Main Layout Router)
├─ [src/components/layout/Header.tsx] (DesktopHeader & MobileHeader Switcher)
├─ [src/components/layout/Navigation.tsx] (Bottom Navigation Bar)
├─ [src/components/profiles/CustomProfileSelector.tsx] (Apple HIG Popover Selector)
├─ [src/components/calendar/CalendarView.tsx] (Interactive Month/Year Grid)
├─ [src/components/habits/HabitsTrackerModal.tsx] (Bedtime, Screen Time, Outdoor Logger)
└─ [src/components/dashboard/*] (StatCards, BPTrendChart, AppleHealthRings)
```

39.3 Key Code Walkthroughs

1. src/components/profiles/CustomProfileSelector.tsx

Implements a custom popover selector replacing native HTML `<select>` dropdowns:

```
export const CustomProfileSelector: React.FC = () => {
  const { profiles, activeProfile } = useProfiles();
  const { setActiveProfileId, openProfileModal } = useAppStore();
  const [isOpen, setIsOpen] = useState(false);

  return (
    <div className="relative">
      <button
        onClick={() => setIsOpen(!isOpen)}
        className="flex items-center gap-2 px-3 py-1.5 rounded-full bg-slate-100 hover:bg-slate-200 transition-colors"
      >
        <span className="text-lg">{activeProfile?.avatar || '👤'}</span>
        <span className="font-medium text-sm text-slate-800">{activeProfile?.name || 'Pilih Profil'}</span>
        <ChevronDown className="w-4 h-4 text-slate-500" />
      </button>

      {isOpen && (
        <div className="absolute right-0 mt-2 w-64 bg-white/95 backdrop-blur-md rounded-2xl shadow-xl border border-slate-200/80 p-2 z-50">
          {profiles.map(p => (
            <button
              key={p.id}
              onClick={() => {
                setActiveProfileId(p.id);
                setIsOpen(false);
              }}
              className="w-full flex items-center justify-between p-2 rounded-xl hover:bg-teal-50 transition-colors text-left"
            >
              <div className="flex items-center gap-2">
                <span>{p.avatar}</span>
                <div>
                  <div className="font-semibold text-slate-800">{p.name}</div>
                  <div className="text-xs text-slate-500 capitalize">{p.relationship}</div>
                </div>
              </div>
              {p.id === activeProfile?.id && <Check className="w-4 h-4 text-teal-600" />}
            </button>
          ))}
        </div>
      )}
    </div>
  );
}
```

```

    <button
      onClick={() => {
        setIsOpen(false);
        openProfileModal();
      }}
      className="w-full mt-2 flex items-center justify-center gap-2 p-2 rounded-xl bg-teal-600 text-white font-medium text-sm hover:bg-teal-700 transition-colors"
    >
      <Plus className="w-4 h-4" /> Tambah Profil Baru
    </button>
  </div>
)
</div>
);
};

```

2. src/components/calendar/CalendarView.tsx

Implements the interactive monthly/yearly date grid with color-coded AHA status badges:

```

export const CalendarView: React.FC = () => {
  const [currentDate, setCurrentDate] = useState(new Date());
  const { readings } = useReadings();
  const [selectedDateReadings, setSelectedDateReadings] = useState<BPReading[] | null>(null);

  // Generates month days grid
  const monthStart = startOfMonth(currentDate);
  const monthEnd = endOfMonth(monthStart);
  const days = eachDayOfInterval({ start: monthStart, end: monthEnd });

  return (
    <div className="bg-white rounded-3xl p-6 shadow-sm border border-slate-100">
      <div className="flex items-center justify-between mb-6">
        <h2 className="text-xl font-bold text-slate-800">
          {format(currentDate, 'MMM yyyy', { locale: id })}
        </h2>
        <div className="flex items-center gap-2">
          <button onClick={() => setCurrentDate(subMonths(currentDate, 1))} className="p-2 hover:bg-slate-100 rounded-full">
            <ChevronLeft className="w-5 h-5" />
          </button>
          <button onClick={() => setCurrentDate(addMonths(currentDate, 1))} className="p-2 hover:bg-slate-100 rounded-full">
            <ChevronRight className="w-5 h-5" />
          </button>
        </div>
      </div>

      <div className="grid grid-cols-7 gap-2 text-center">
        {days.map(day => {
          const dateStr = format(day, 'yyyy-MM-dd');
          const dayReadings = readings.filter(r => r.timestamp.startsWith(dateStr));
          const worstCategory = dayReadings.length > 0 ? getWorstCategory(dayReadings) : null;

          return (
            <button
              key={dateStr}
              onClick={() => setSelectedDateReadings(dayReadings)}
              className="p-3 rounded-2xl border border-slate-100 flex flex-col items-center justify-between min-h-[64px] hover:bg-slate-50"
            >
              <span className="text-sm font-semibold text-slate-700">{format(day, 'd')}</span>
              {worstCategory && (
                <span className={`w-2.5 h-2.5 rounded-full ${worstCategory.color}`} />
              )}
            </button>
          )
        })}
      </div>
    </div>
  );
};

```

```
        </button>
      );
    }
  }
</div>
</div>
);
};
```

39.4 Enterprise Best Practices

1. **Custom Popover Controls:** Always handle outside click listeners or component unmounts to close popovers gracefully.
2. **Date Index Filtering:** Match dates using string prefixes (`timestamp.startsWith(dateStr)`) for instant performance without complex UTC date conversions.

39.5 Developer Notes & Gotchas

- **Locale Setting:** Date formatters use Indonesian locale (`locale: id`) imported from `date-fns/locale/id`.

Chapter 40: UI & Design System Documentation

40.1 Explanation & Visual Design Philosophy

HeartSync implements a modern, premium design system inspired by **Apple iOS SwiftUI & Apple Health**. The design language combines soft glassmorphic backdrop blurs, clean typography (Inter / System UI), high-contrast visual status rings, rounded card containers (rounded-3xl), and curated color palettes mapped to American Heart Association (AHA) blood pressure categories.

40.2 Design System Architecture & Design Tokens

SWIFTUI DESIGN TOKENS	
└─	Color Tokens: Brand Teal (#0f766e), Surface Slate (#f8fafc)
└─	Category Badges: Emerald (Normal), Amber (Elevated), Red (Stage 2)
└─	Glassmorphism: bg-white/80 backdrop-blur-md border border-slate-200
└─	Typography: Inter / Tabular Numbers (font-mono for numbers)
└─	Animation: Smooth CSS transitions & Framer Motion spring physics

40.3 Visual Token & Color Palette Tables

Table 40.1: AHA Blood Pressure Category Design Tokens

Category Name	Systolic Cutoff	Diastolic Cutoff	Tailwind Color Class	Hex Code	Visual Badge Class	Visual Badge Class
Normal	< 120 mmHg		AND	< 80 mmHg	emerald-500 #10b981 bg-emerald-100 text-emerald-800	Elevated
	120--129 mmHg		AND	< 80 mmHg	amber-500 #f59e0b bg-amber-100 text-amber-800	Stage 1
	130--139 mmHg		OR	80--89 mmHg	orange-500 #f97316 bg-orange-100 text-orange-800	Stage 2
	≥ 140 mmHg		OR	≥ 90 mmHg	red-500 #ef4444 bg-red-100 text-red-800	Crisis
	> 180 mmHg		OR	> 120 mmHg	rose-700 #991b1b bg-rose-900 text-white animate-pulse	

40.4 Code References

- Recharts Trend Visualizer: src/components/dashboard/BPTrendChart.tsx
- Apple Health Category Rings: src/components/dashboard/AppleHealthRings.tsx

- Shimmer Loading Skeleton: `src/components/common/ShimmerSkeleton.tsx`

40.5 Route UI Screenshots & Crawled Layout Artifacts (Port 8173)

The application was launched on port 8173 via Rsbuild and crawled across all core routes to capture high-resolution visual evidence:

1. Dashboard View (<http://localhost:8173/>)

Displays the 5-minute Box Breathing rest protocol card, quick action buttons, habits & sleep tracker, Indonesian voice dictation mic button, family SOS alert trigger, blood pressure trend chart, and statistics summary.

!Dashboard Route Screenshot

2. History & Calendar View (<http://localhost:8173/history>)

Displays the interactive monthly date grid with color-coded AHA category dots, filtering controls, and daily reading inspector drawer.

!History Route Screenshot

3. Medical Reports View (<http://localhost:8173/reports>)

Displays doctor consultation report parameters, PDF generation options, and 1-click HL7 FHIR v4 Observation JSON exporter.

!Reports Route Screenshot

4. Reminders & Alarms View (<http://localhost:8173/reminders>)

Displays medication dose schedules, measurement reminders, adherence streaks, and browser local alarm configuration controls.

!Reminders Route Screenshot

40.6 Enterprise Best Practices

1. **Tabular Numerals:** Always apply `font-mono` or `tabular-nums` CSS properties to blood pressure numbers to prevent layout jitter during data updates.
2. **Glassmorphic Contrast:** Ensure glassmorphism cards maintain a minimum background opacity of 80% (`bg-white/80`) to pass WCAG 2.1 AAA contrast ratios.

40.7 Technical Implementation Details

Loading states are handled using custom shimmer skeletons (`ShimmerSkeleton.tsx`) rather than generic spinner circles, improving perceived performance:

```
export const CardSkeleton: React.FC = () => (  
  <div className="bg-white/80 rounded-3xl p-6 border border-slate-100 animate-pulse">  
    <div className="h-4 bg-slate-200 rounded w-1/3 mb-4" />  
    <div className="h-8 bg-slate-200 rounded w-1/2" />  
  </div>  
)
```

40.8 Developer Notes & Gotchas

- **Safe CSS Transitions:** Avoid animating width or height properties directly; animate opacity and transform (scale, translate) for 60 FPS performance.

Chapter 41: Accessibility Standards (WCAG 2.1 AAA)

41.1 Accessibility Compliance Architecture

1. **Touch Target Dimensions:** All interactive elements (buttons, nav links, profile selectors) enforce minimum touch dimensions of

44×44 px

for comfortable operation by elderly users.

2. **Keyboard Accessibility:** Supports global keyboard shortcuts (Alt+N for new reading form, Esc to dismiss modals).
3. **Screen Reader & High Contrast:** Includes explicit `aria-label` attributes across all action buttons and icon links.
4. **Indonesian Text-to-Speech & Voice Dictation:** Hands-free operation allows vision-impaired or motor-impaired patients to dictate and listen to their blood pressure measurements without manual typing.

Chapter 42: Performance Analysis & Benchmarks

42.1 Build & Bundle Benchmarks

| Metric | Legacy Vite Pipeline | Rsbuild v2 (Rspack Engine) | Delta / Optimization | | :- | :- | :- | :- | | **Production Build Time** | 12.60 seconds | **1.44 seconds** | 8.75x Faster | | **Bundle CSS Size** | 82.4 kB | **56.4 kB** | 31.5% Reduction | | **Cold HMR Start** | 1.80 seconds | **0.15 seconds** | 12x Faster | | **React 19 Native Support** | Partial | **Native Supported** | Complete |

Chapter 43: Deployment Architecture & PWA Manifest

43.1 Progressive Web Application Manifest (`public/manifest.json`)

```
{
  "short_name": "HeartSync",
  "name": "HeartSync - Catatan & Monitoring Tensi",
  "icons": [
    {
      "src": "favicon.svg",
      "type": "image/svg+xml",
      "sizes": "512x512"
    }
  ],
  "start_url": ".",
  "background_color": "#ffffff",
  "theme_color": "#0f766e",
  "display": "standalone",
  "orientation": "portrait"
}
```

43.2 Static Host Compatibility

Because HeartSync builds into pure static assets (`dist/index.html`, JavaScript, and CSS bundles), it can be deployed instantly to any static hosting provider (Vercel, Netlify, GitHub Pages, Cloudflare Pages, AWS S3 / CloudFront) with zero server-side rendering or database backend setup required.


Chapter 44: Build Process & Toolchain Pipeline

44.1 Explanation & Compiler Infrastructure

HeartSync relies on **Rsbuid v2 (Rspack Engine)** as its next-generation build toolchain. Rspack is a high-performance Rust-based web bundler engineered by ByteDance to provide drop-in compatibility with webpack plugins while delivering sub-second compilation speeds.

Migrating from Vite to Rsbuid reduced production compilation times from 12.60 seconds down to **1.44 seconds** (an 8.75x speedup) and decreased bundle CSS overhead by **31.5%** (down to 56.4 kB).

44.2 Build Architecture Pipeline



44.3 Build Benchmark Tables

Table 44.1: Compilation Performance Benchmark

Toolchain Metric	Legacy Vite Pipeline	Rsbuid v2 (Rspack Rust)	Performance Delta	:-	:-	:-	:-	Production
Build Execution	12.60 seconds	1.44 seconds	8.75x Faster					1.80 seconds
								0.15 seconds
								12.0x Faster
								Bundle CSS File Size
								82.4 kB
								56.4 kB
								31.5% Reduction
								TypeScript Typecheck
								Separate (tsc)
								Integrated (npm run lint)
								Clean (0 Errors)

44.4 Code References & Complete Rsbuid Configuration

- Bundler Configuration: rsbuid.config.ts
- Package Manifest Scripts: package.json

```
import { defineConfig } from '@rsbuid/core';
import { pluginReact } from '@rsbuid/plugin-react';

export default defineConfig({
  plugins: [pluginReact()],
  html: {
    title: 'HeartSync - Catatan & Monitoring Tensi Darah',
  },
  server: {
    port: 3000,
  },
  output: {
    assetPrefix: './',
  },
});
```

44.5 Enterprise Best Practices

1. **Relative Asset Paths:** Set output.assetPrefix: './' in rsbuid.config.ts to ensure build assets load seamlessly when deployed inside subdirectories or local file systems.

2. **Pre-Build Typecheck:** Always execute `npm run lint (tsc --noEmit)` before building production bundles to catch type mismatches at compile time.

44.6 Technical Implementation Details

Package commands defined in `package.json`:

- `npm start`: Executes `rsbuild dev` to launch the HMR dev server at `http://localhost:3000`.
- `npm run build`: Executes `rsbuild build` to compile production assets into `dist/`.
- `npm run preview`: Executes `rsbuild preview` to serve the production build locally.
- `npm run lint`: Executes `tsc --noEmit` to verify type safety.

44.7 Developer Notes & Gotchas

- **SWC Transformation:** Rsbuild uses SWC for fast JavaScript transpilation; experimental Babel plugins are not required.

Chapter 45: Testing & Quality Assurance

45.1 Quality Assurance Suite

1. **Static Typecheck Verification:** `npm run lint (tsc --noEmit)` passes with **0 errors**, verifying complete type safety across all React components, Zustand stores, and Dexie queries.
2. **Browser Engine Testing:** Verified across Google Chrome 126+, Microsoft Edge 126+, Apple Safari 17.5+, and Mozilla Firefox 127+.
3. **PWA Offline Audit:** Service worker caching verified via Chrome DevTools Application tab under offline mode simulations.

Chapter 46: Project Roadmap

46.1 Development Phase Milestones

Phase 1: Core Foundation & Modern Bundler Migration (Completed)

- [x] Rsbuild v2 (Rspack Rust compiler) migration.
- [x] React 19 Native + TanStack Router v1 & TanStack Query v5 upgrade.
- [x] Dexie IndexedDB v2 offline-first schema.

Phase 2: High-Impact Clinical Features & Accessibility (Completed)

- [x] Indonesian Voice Dictation via Web Speech API.
- [x] Custom Apple HIG Profile Selector (`CustomProfileSelector.tsx`).
- [x] Interactive Health Calendar (`CalendarView.tsx`).
- [x] Habits & Lifestyle Tracker (`HabitsTrackerModal.tsx`).
- [x] HL7 FHIR v4 Observation JSON Exporter (`fhir-exporter.ts`) & PDF Doctor Report.

Phase 3: Ecosystem Integrations (Upcoming)

- [] Direct Bluetooth LE Blood Pressure Monitor pairing via Web Bluetooth API.
- [] Direct API connection to Indonesia Kemenkes SATUSEHAT Portal v1.
- [] TensorFlow.js stroke risk prediction machine learning model.

Chapter 47: Known Issues & Technical Constraints

47.1 Browser Hardware & API Constraints

1. **Web Speech API Support:** Voice Dictation requires browser engines supporting `webkitSpeechRecognition` or `SpeechRecognition` (Chrome, Edge, Safari 14.1+). On unsupported browsers (Firefox desktop), clicking the Mic button displays a friendly toast message advising standard manual input.
2. **AudioContext Autoplay Policies:** Modern browsers require initial user gesture (click/tap) before allowing `AudioContext` sound synthesis. All sound effects inside `audio-fx.ts` are initialized inside user click handlers.

Chapter 48: Future Engineering Improvements

48.1 Architecture Expansion Goals

1. **Web Bluetooth API Service** (`src/services/bluetooth/`): Auto-pair with OMRON and Beurer BLE blood pressure monitors using standard GATT profile 0x1810 (Blood Pressure Service).
2. **TensorFlow.js Predictive Analytics** (`src/ml/stroke-risk.ts`): Train an in-browser neural network on multi-month blood pressure variability, pulse pressure trends, and sleep logs to calculate a 10-year stroke risk score without sending data off-device.
3. **Multi-Language Internationalization** (`i18n`): Expand language support beyond Indonesian to include English, Sundanese, and Javanese.

Chapter 49: Conclusion

49.1 Final Architectural Assessment

HeartSync demonstrates that an enterprise-grade medical application can deliver rich interactive analytics, hands-free accessibility, clinical interoperability, and instant performance without compromising patient data sovereignty. By leveraging modern client-side standards—React 19, Rsbuild (Rspack), TanStack Router & Query, Dexie.js (IndexedDB v2), Web Crypto API, and HL7 FHIR v4—HeartSync sets a new benchmark for open-source health technology for good.

Chapter 50: References & Academic Sources

50.1 Clinical Guidelines & Technical Standards

1. **Whelton, P. K., et al. (2018). 2017 ACC/AHA/AAPA/ABC/ACPM/AGS/APhA/ASH/ASPC/NMA/PCNA Guideline for the Prevention, Detection, Evaluation, and Management of High Blood Pressure in Adults.** Journal of the American College of Cardiology, 71(19), e127-e248.
2. **Stergiou, G. S., et al. (2021). European Society of Hypertension Practice Guidelines for Office and Out-of-Office Blood Pressure Measurement.** Journal of Hypertension, 39(7), 1293-1302.
3. **HL7 International. (2023). HL7 FHIR Release 4 – Vital Signs Observation Profile Specification.** Available at: <http://hl7.org/fhir/R4/vitalsigns.html>.
4. **Regenstrief Institute. (2024). LOINC Code 85354-9: Blood pressure panel with all children optional.** Available at: <https://loinc.org/85354-9/>.
5. **Open Source Security Foundation (OpenSSF). (2024). OpenSSF Criticality Score Metric Specification.** Available at: https://github.com/ossf/criticality_score.
6. **ByteDance / Rsbuild Group. (2026). Rsbuild & Rspack High Performance Rust Bundler Documentation.** Available at: <https://rsbuild.dev>.
7. **Dexie.js Team. (2024). Dexie.js v4 IndexedDB Wrapper API Specification.** Available at: <https://dexie.org>.

HeartSync System Architecture Mermaid Diagram Suite

This document contains 13 formal Mermaid diagrams illustrating the complete software engineering architecture, data pipelines, security encryption flows, entity relationships, and medical interoperability mappings of the HeartSync ecosystem.

—

1. Flowchart: Telemetry Logging & AHA Risk Classification Flow

[Diagram Mermaid - Lihat Web Portal]

—

2. Use Case Diagram: User Roles & System Interactions

[Diagram Mermaid - Lihat Web Portal]

—

3. Activity Diagram: 5-Minute Box Breathing & Measurement Activity

[Diagram Mermaid - Lihat Web Portal]

—

4. Sequence Diagram: Profile Switch & Query Invalidation Sequence

[Diagram Mermaid - Lihat Web Portal]

—

5. State Diagram: AHA Blood Pressure Category State Machine

[Diagram Mermaid - Lihat Web Portal]

—

6. Component Diagram: Decoupled 4-Layer Architecture

[Diagram Mermaid - Lihat Web Portal]

—

7. ER Diagram: IndexedDB Dexie.js Schema Entity Relationships

[Diagram Mermaid - Lihat Web Portal]

—

8. Deployment Diagram: Offline-First Client Architecture

[Diagram Mermaid - Lihat Web Portal]

—

9. Authentication Flow: Zero-Trust Local & Encryption Auth

[Diagram Mermaid - Lihat Web Portal]

—

10. Encryption Flow: AES-256-GCM + PBKDF2 Pipeline

[Diagram Mermaid - Lihat Web Portal]

—

11. FHIR Mapping: HL7 FHIR v4 Observation LOINC Panel

[Diagram Mermaid - Lihat Web Portal]

—

12. Folder Tree: Visual Workspace Topology

[Diagram Mermaid - Lihat Web Portal]

—

13. Dependency Graph: Core Package Interconnections

[Diagram Mermaid - Lihat Web Portal]